



## Quick Coomassie Stain User Guide

Quick Coomassie is a new revolution in rapid 1-step Coomassie staining. The proprietary formulation, incorporating Colloidal Coomassie, is used for rapid protein staining in polyacrylamide gels.

### Storage Conditions:

Upon receipt, store the stain at Room Temperature. The QC stain is stable for a minimum 12 months at Room Temperature.

### Simple 1-step Protocol:

1. Pour 25 ml QC stain into a container. Use more stain if you are using a larger gel tray.
2. Remove the gel from the cassette and place the gel into the stain.
3. Leave the gel, while shaking, for a minimum of 15 minutes or until all weak protein bands are fully developed. Stain intensity is high after about 1 - 2 hours and maximum after overnight incubation.
4. Transfer the gel to DI water to remove any background staining and for gel storage. (N.B: A minimum 1 hour full stain is recommended before storing the gel in water.)

### Microwave Procedure for Gels:

1. Using a microwave to heat up the QC stain can speed up the development of the protein bands.
2. For turbo-charging the stain, we recommend microwaving the gel, immersed in QC stain, in a suitable microwave-safe tray for a maximum 10 seconds at full power.
3. Remove the tray from the microwave and keep the gel in the QC stain for at least 30 min – 1 hour before storing the gel in DI water.

### For Mass Spectrometry Applications:

1. Stain the gel as normal.
2. Excise the protein band of interest and put in a clean microfuge tube ideally.
3. Add 1 ml of 30% ethanol or 30% acetone.
4. Incubate for 20 min (60°C – 70°C increases the rate of de-staining).
5. Decant supernatant and repeat step 3 and 4 at least 3 times or until the gel fragment is clear.
6. Run the sample in the mass spectrometer.

### Ordering Information

Product	Units	Order Code
1 Litre QC Stain (sufficient for 40 gels if used once)	1 L	GEN-QC-STAIN-1L
3 Litre QC Stain (sufficient for 120 gels if used once)	3 x 1 L	GEN-QC-STAIN-3L

### Protein Ark Limited

Telephone: +44 (0) 33 33 44 20 25

FAX: +44 (0) 33 33 44 20 25

Email: [info@proteinark.com](mailto:info@proteinark.com)

### In association with Generon Limited

**Generon**

